

Hormone Profiling: Saliva, Blood & Urine

Dr Aron Gonshor: PhD, DDS, FRCD(C), FAO

Accurately Assessing Sex and Adrenal Hormones

- Which 'Medium' should be used?
- Which is the 'Best'?
 - What questions need answers?

Consideration of the Relative Merits of the Three Primary Analytical Mediums Used to Test Hormones

Blood (Serum)

Blood

Bound + Free Steroids

(~97%)

(~3%)

Albumin
(Non-specific)

SHBG

(Sex Hormone Binding Globulin)
Binds Estrogens & Androgens

CBG

(Cortisol Binding Globulin)
Binds Cortisol & Pg

Steroid Hormones in Blood

Serum Sampling

Advantages

- Readily Available
- Cost-effective
- Widely utilized and accepted
- Well-established reference ranges
- Measures SHBG & CBG
- Suitable for testing peptide hormones
 - eg., Prolactin, TSH, FSH, LH:
Due to their molecular weight, do **not** appear in filtrates, in saliva or urine
 - Mostly reimbursable by insurance

Disadvantages

- **Spot Test:** Provides snapshot in time
 - Hormone levels time-dependent, with ultradian, circadian & diurnal periodicity
- **Poor reliability for pulsatile secretion: Major fluctuations, especially in sex hormones**
- Short half-life. Rapid clearance & compartmentalization
- **Generally no distinction between bound & free hormone fractions**
- **Invasive collection**
- Specificity often poor in immunologically based assays; eg., RIA & ELISA
- Full hormone profiling is cost-prohibitive.


SHBG = Sex Hormone Binding Globulin; CBG = Cortisol Binding Globulin

Ultradian = recurrent cycles in a 24-hour day; Circadian = Complete one cycle daily

Consideration of the Relative Merits of the Three Primary Analytical Mediums Used to Test Hormones

Saliva

Bound vs Free Steroid Fractions



Saliva Sampling

Advantages

- Excellent for ***Rapid testing***
 - Accurate measure of validated biomarkers
- Well-established reference ranges
- **Non-invasive**, with low infection rate
- Measures **free hormone fraction**
- Ideal for **multiple or serial sampling**
 - More versatility than serum for evaluating un-supplemented hormone status
- Accessible to practitioners, without a license, to order blood tests

Disadvantages

- **Production** difficult for some pts
 - Eating, drinking, gum-chewing, use of make-up
 - Inadequate volume in some > 40 yo (> females)
 - (Sjogren's Syndrome, etc)
 - Micro-damage from tooth brushing (~1hr prior to testing), even with no bleeding signs
- **No measure of hormone metabolites**
- Collection can, itself, significantly alter hormone concentrations

Salivary Hormone Testing

- In bloodstream: Endogenously or exogenously-produced hormones
 - 97–98% of the active hormone are bound up by specific proteins
 - **2–3%** released in capillary beds and into interstitial space & tissues
- In saliva - most metabolites from oral hormone therapy filtered out by salivary gland. Only **active** hormones enter saliva
 - ***Salivary hormones more representative of the amount of **bioactive steroid** present in the bloodstream, & its **bioavailability** to target tissues***

FLUIDS iQ - Stress and Hormone Profile

17- β ESTRADIOL (E2) ng/ml

4.9

CORTISOL (C) ng/ml

Stress and Hormone Profile: Standard & Advanced

A) Standard:

- **Adrenal function, Neuropsychiatric Disorders & Aging:**
 - Diurnal Cortisol, DHEA
- **Sex Hormones:** Estrogen, Progesterone, Testosterone
- **Intestinal Immunity:** Secretory Immunoglobulin (SIgA)
- **Circadian Rhythm & Anti-Oxidant Status:** Melatonin

B) Advanced: Addition of the following

- **Metabolic Function:** Insulin – Fasting & Non-Fasting

Mid-Range

120 - 220

High

223 - 359

Very High

> 359


Greater than 69

Range not applicable

Consideration of the Relative Merits of the Three Primary Analytical Mediums Used to Test Hormones

Urine

Diurnal Pattern of Testosterone (T) Secretion




Brenner et al, 1983

- Testosterone levels:
 - Peak at ~ 8 AM (0800), then **fluctuate constantly** throughout the day
 - Lowest at ~ 8 PM (2000), then climb at night
- In 40yr old, AM readings may be 200 points higher than in evening, versus a 50-point difference for a 70 year old

- Interpretation of results of hormone measurements requires recognition of the biological daytime and day-to-day changes of hormones. [Ahokoski, 1998](#)

Daytime Variation in Serum Progesterone (Pg) During Mid-Luteal Phase



- Pg released into bloodstream in **pulsations**, leading to large fluctuations in serum Pg.
- Up to **8** fold in 90 min & **2.3- 40.1ng/ml (8-140 nmol/L)**, or ~ **20** fold change over 24 hrs

Single blood or urine test unreliable

Ideally, test blood or urine over 24 hrs, or 2-3 days in a row, to get average*

*Beckley, 2021

24 *hr* Urine Sampling

• Advantages:

- Shows **unbound** (bio-available) fraction
- Allows determination of **metabolites**
- Prevents **false** highs & lows as occurs when testing '**spot**' methods like serum, as well as dried urine
- Clinically:
Correlates well with **symptoms** on hormone questionnaires
 - Allows evaluation of adequacy & safety of **exogenous estrogen supplementation**

• Disadvantages:

- Results affected by significant liver or renal impairment, dehydration or excessive fluid intake: **Fundamental Creatinine issue**
- Result misinterpretation possible due to hepatic '1st-pass' effect, if on oral HRT
- **Does not elucidate *diurnal cortisol***
 - Best in salivary testing

Urine Testing & Hormone Supplementation

- Urine testing reveals the ~ **total hormone consumed & eliminated**, but **NOT** how much of the **bioactive & available** hormone is present in the systemic circulation, or has already entered target tissues.
- **Conclusion:** Urine testing an excellent way to evaluate **endogenous** production of the sex hormones, but:

Not clinically useful as a diagnostic fluid for exogenous oral or topical hormone delivery

***Dried Urine* Spot Sampling**

- All samples obtained at a defined moment in time are, by definition:
A SPOT TEST
- A Spot test is a '**snapshot in time**', which can be useful for certain analytes, but has fundamental flaws & compromises, especially:
 - - for pulsatile analytes (eg; hormones)
 - - when compared to a full 24-hr hormonal production cycle
- Forced to normalize to the breakdown product of creatine phosphate from muscle and protein metabolism, & found in urine - ***Creatinine***

Dried Urine & Creatinine Normalization

Ideally, want to know how much hormone is produced in the classic $\mu\text{g}/24 \text{ hrs}$.
This eliminates urine volume variability (within day & day-to-day)

- Instead, forced to normalize to **Creatinine (Cr)**, shown as ng of analyte/mg Cr
 - Fundamental problem for **pulsatile** analytes, such as hormones
 - Creatinine levels vary widely with age, sex, diet, exercise, kidney function & genetics
 - Subject ends up as own control. **Cannot be classified with a typical reference range**
- Reference ranges quoted for spot urine samples span an order of magnitude, that is a **factor of 10!**
- In contrast, span for 24-hr urine sample lies between a factor of **2 to 3**

Dried Urine (DU) Spot Tests

- In absence of independently published methods/data & given the limited urine volume in a Dried Urine sample for this type of analysis, it is unclear how quantitation of certain analytes can be **credibly accomplished**, **accurately assessed** and **reproduced** with confidence

The Dried Urine paradigm is arguably outside acceptable scientific standards

Serum Testing for Hormones

- Limited utility for sex hormones, because:
 - No distinction made between **bound & free** hormone.
 - Estradiol, estrone, estriol and Pg reported as **total** hormones
 - **Free** hormone assays not common, leading to misleading results
 - Hormone levels appear to be normal, or even high normal, because of an abundance of **bound** hormone

However, if the **free** hormone level is **low**, patient can be **functionally deficient**, even with a **normal total** hormone level

Appropriate Body Fluids for Testing **Exogenously** Delivered Hormones

Type of Fluid	Non Endogenous	Oral Steroids	Topical Steroids	Vaginal Steroids	Troche Steroids	Pellet/IM Steroids
Serum	Yes	Yes ¹	No ²	No ²	Yes	Yes
Saliva	Yes	Yes	Yes ³	Yes	No ⁴	Yes
Urine	Yes	Yes ¹	No ²	No ⁵	Yes	Yes
DBS	Yes	Yes ¹	Yes ⁶	Yes	Yes	Yes

1. Overestimation: Possible metabolite interference with immunoassays (GC/MS OK)

Saliva & DBS, but not serum/plasma or urine, are the only way to **accurately** monitor exogenous **Topical** and **Vaginal** steroid hormone therapies

6. Overestimation: If hands/fingertips used to apply hormones within 24 hrs of testing

Conclusions: Fluid Mediums For Hormone Testing

For Exogenous Hormones

- **Serum:** For accurate oral administration of sex hormones, needs extraction & separation of metabolites
- **Urine:** Not recommended unless range re-established for dosing (eg, 100–300 mg oral Pg) & range readjusted up to reflect the expected level
- **Serum/plasma & Urine:** NOT recommended for testing steroid hormones delivered by **topical** route of administration.

Best with Saliva

Hormone Profiling: Saliva, Blood & Urine

Dr Aron Gonshor: PhD, DDS, FRCD(C), FAO